


— 11.458

8.184
7.591
7.571
7.403
7.390
7.370
7.230
5.864
4.307
4.295
3.849
3.823
3.592
3.582
3.571
3.491
3.322
3.285
3.257
3.229
3.104
3.087
3.069
2.508
2.372
2.255
2.215
2.113
1.685
1.657



NAME EPZ-6438
EXPNO 1
PROCNO 1
Date_ 20170915
Time 15.42
INSTRUM spect
PROBHD 5 mm PABBO BB-
PULPROG zg30
TD 65536
SOLVENT DMSO
NS 8
DS 0
SWH 8223.685 Hz
FIDRES 0.125483 Hz
AQ 3.9846387 sec
RG 144
DW 60.800 usec
DE 6.50 usec
TE 300.0 K
D1 1.00000000 sec
TD0 1

==== CHANNEL f1 =====
NUC1 1H
P1 10.85 usec
PL1 -4.00 dB
PL1W 20.19063568 W
SFO1 400.1324710 MHz
SI 32768
SF 400.1300000 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00