




Sample Information

Sample Name:	3X FLAG Peptide	Gas Temp:	250 °C
Operator:	Jefferson	Gas Flow:	10L/min
Injection Date:	12/23/2017 4:07:12 PM	Nebulizer:	50psi
Batch No.:	GP101491	VCap:	4500v
T.Flow:	0.2ml/min	Fragmentor:	70v

Result



Molecular Weight:	2861.87
m/z([M+3H]3+) Base Peak:	954.90
m/z([M+4H]4+) Base Peak:	716.40
m/z([M+5H]5+) Base Peak:	573.30
m/z([M+6H]6+) Base Peak:	477.85